



ENTA[®]
Ceramic

Vormenkaart


Mould chart
Formenkarte
Carte de formes

ceramic

Upper anteriors



Upper anteriors



Upper anteriors



Lower anteriors



Lower anteriors



ceramic


Posteriors



- All moulds are illustrated actual size.
- Available in Ao-D4, including Bo
- The dimensions are in millimetres.

Optiform posteriors (Lingualized occlusion concept)

ceramic



combination table ENTA® assortment

	ANTERIOR TEETH			POSTERIOR TEETH			
	Upper	Lower E/H	Lower LOC	20° (E)	10° (H)	Optiform *	
Oval	110	62 - 70	70 - 72	E2	H4	34S - 30	
	120	62	70 - 72	E4 - E5		34	
	130	70	70 - 72	E2	H4	34S - 30	
	140	72	76	E4		34	
	150	76	78	E2 - E4	H4	34S - 30	
Rectangular	210	64	72	E4		34	
	220	64	72	E4		34	
	230	74	76	E5		34	
	240	66 - 80	80 - 84	E4		34	
	250	76	84	E4 - E5		34	
	260	82	84 - 88	E5		34	
	270	82	84 - 88	E5		34	
	280	86	86 - 88	E5		34	
	Triangular	310	60	62	E2	H4	34S
		320	68	76	E2	H4	34S
330		64	72	E4		34	
340		74	76	E4 - E5		34	
350		70	70 - 72	E2	H4	34S - 30	
360		64 - 68	72	E2 - E4	H4	34S - 30	
370		68 - 76	76	E4 - E5		34	
380		72 - 76	76 - 84	E4		34	
Square	410	62	70 - 72	E4		34 - 34S	
	420	74	76 - 84	E4		34 - 34S	
	430	64	72	E4		34 - 34S	
	440	80	82	E5		34	
	450	76 - 82	84	E5		34	
	460	74	76 - 84	E4		34 - 34S	

* Important, see also manual L.O.C.